LADY <u>Reques</u>	Lab reference number			
Name	:	BHT/Clinic No :		
Age	:	Ward/Clinic :	•••••	
Gender	:M 🗌 F 🗌	Tel No of patient :		
Date of birth	:	Specimen : VENOUS IN 2CC E	BLOOD DTA TUBE	
Fasting status : <4 Hrs  4 Hrs  6 Hrs				

	 	 	 •••
Drugs :	 	 	 •••

## PLASMA



## Homocysteine

Name and designation of the requesting medical officer :	To be filled by the nursing officer Date & time of sample collection :	
Date :	Name :	
To be filled at the laboratory reception   Date & time of sample receipt   Lipae- mic   Hae- mic	molyzed 🔲 Icteric 🔲 Insufficient	

Please follow the instructions provided overleaf. For further information and clarifications please contact Department of Chemical Pathology 011–2693711 Ext 314

## Instructions for patient preparation and sample collection

## PLASMA HOMOCYSTEINE

Patient Preparation - Fasting	
Children > 3 years	: 4-6 hrs
Children $> 1$ yr $- 3$ years	: 3 hrs
Infants	: Just before the next feed
මේම් පරික්ෂණය සඳහා රෝගියා	මරාතාරව සිටිය යුතු කාල සිමාව
අවුරුදු 3 ට වැඩි ළමුන්	: <b>5763</b> 4 - 6
අවුරුදු 1 ත් 3 ත් අතර ළමුන්	: 5763 3
අවුරුද්දට අඩු	: ඊළහ ආහාර/කිරි වෙිලට පෙර
ஆகாரமின்றி இருப்பதற்கான த	யார்ப்படுத்தல்
சிறுவர் >3 வயது	: 4-6 hrs
சிறுவர் > 1 வயது - < 3 வயது	: 3 hrs
குழந்தை	் அடுத்த தாய்ப்பால் கொடுக்க சற்று முன்

Collection Container	Tube containing EDTA	
Sample Volume	2 mL	
Collection and transport instruc- tions	Sample should be taken without a tourniquet. Once the sample is transferred to the tube, mix it well and transport it in the ice to the laboratory immediately.	